



# Capillary Automated Footprinting Analysis



Somdeb Mitra, Inna Shcherbakova, Russ B. Altman, Michael Brenowitz and Alain Laederach  
<http://simtk.org/cafa>. E-mail: Alain Laederach (alain@helix.stanford.edu)

User's Manual for the CAFA software.  
March 2008, Stanford University  
Alain Laederach

## Installation

Platform specific installation instructions are provided in separate documents that can be downloaded from <https://simtk.org/home/cafa>. The figures in this manual are from the Mac Intel version of CAFA, but are equivalent to the Windows version.

## General Data Philosophy

CAFA is currently designed to read in output from a Beckman CEQ-8000 sequencer. The data has to be exported from the Beckman software using the "Text Export" functionality in the database. The output files should look like this:

Sample Name: 7.A03\_07032003TN  
Sample Subject ID:

### Raw Data Output Injection

INDEX	CAP	FILTER 1	FILTER 2	FILTER 3	FILTER 4	CURRENT	VOLTAGE	RAW	CURR				
1	A	636	346	804	242	1.39	0.0	0.00	0.00	0.00	0.0	0.00	0
2	A	634	356	764	236	1.52	1.1	0.00	0.00	0.00	1.1	0.00	0
3	A	558	318	816	242	1.64	0.9	7.00	0.00	0.00	0.9	0.00	0
4	A	564	320	788	196	1.75	1.8	11.00	23.00	2.50	1.8	23.00	18
5	A	590	314	812	250	1.86	1.9	9.00	23.00	3.38	1.9	23.00	18
6	A	592	324	844	222	1.97	1.9	8.00	27.00	3.38	1.9	27.00	26
7	A	562	360	746	256	2.06	1.9	7.00	27.00	3.00	1.9	27.00	26
8	A	626	304	840	248	2.14	1.9	10.00	26.00	2.88	1.9	26.00	26
9	A	586	336	744	260	2.21	1.9	9.00	26.00	3.63	1.9	26.00	26
10	A	574	344	780	262	2.27	1.9	7.00	26.00	3.13	1.9	26.00	26
11	A	648	334	780	256	2.32	1.9	10.00	26.00	2.63	1.9	26.00	26
12	A	626	310	724	246	2.35	1.9	9.00	26.00	3.75	1.9	26.00	26
13	A	598	336	814	224	2.37	1.9	8.00	27.00	3.13	1.9	27.00	26
14	A	600	332	750	232	2.38	1.9	10.00	27.00	2.75	1.9	27.00	26
15	A	616	346	784	242	2.37	1.9	9.00	27.00	3.50	1.9	27.00	26
16	A	594	304	828	216	2.34	1.9	8.00	27.00	3.38	1.9	27.00	26
17	A	604	302	782	238	2.30	1.9	9.00	27.00	3.00	1.9	27.00	26
18	A	608	338	808	250	2.25	1.9	9.00	27.00	3.25	1.9	27.00	26
19	A	580	318	812	244	2.19	1.9	8.00	27.00	3.00	1.9	27.00	26
20	A	574	316	744	234	2.11	1.9	10.00	27.00	3.00	1.9	27.00	26
21	A	584	332	750	246	2.03	1.9	6.00	27.00	3.63	1.9	27.00	26
22	A	628	336	788	224	1.93	1.7	3.00	27.00	3.38	1.7	27.00	26
23	A	620	316	828	208	1.82	0.0	0.00	27.00	4.38	0.0	27.00	26
24	A	554	322	794	230	1.71	0.0	1.00	0.00	0.00	0.0	0.00	0
25	A	562	300	788	262	1.59	0.0	0.00	0.00	0.00	0.0	0.00	0
26	A	630	324	772	240	1.47	0.0	0.00	0.00	0.00	0.0	0.00	0
27	A	620	362	756	242	1.35	0.0	0.00	0.00	0.00	0.0	0.00	0
28	A	632	356	862	256	1.22	0.0	0.00	0.00	0.00	0.0	0.00	0

Currently, we do not directly support other types of input, but depending on need, please contact Alain Laederach ([alain@helix.stanford.edu](mailto:alain@helix.stanford.edu)) if you would like to see other input formats supported.


The basic idea behind the CAFA data philosophy is that each well from the 96 well plate has one corresponding text file that contains the raw trace (Filter 1) and the corresponding reference ladder (Filter 2). It is important when setting up

experiment to always remember to add a reference ladder to the sample. Without the reference ladder, CAFA will not be able to analyze your data and will crash.

As such, when setting up to analyze data, you will want to create a directory with all your data in it. You may want to setup a directory structure for each experiment to avoid an overwhelming amount of data being analyzed simultaneously. For example, if you are mapping the structure of an RNA under various solution conditions, or performing a titration, you may want to have a directory with all the data pertaining to a particular titration in /Titration1, and another Titration in /Titration2. CAFA will generate an output file with all the data in it in the end, and you can organize what data goes in that output file by organizing your data into different directories prior to analysis.


### Reading in the Data

For this example we will consider a directory that has the following three traces (DMS maps of the L-21 group I intron), this data is also available for download from <http://simtk.org/home/cafa>.




The first two files are repeats of the experiment and the seq1\_070711215M.txt file is control lane we will use as background for identifying RT stops.

Open CAFA (launching instructions differ depending on whether you are using a Mac or a PC) and select the working directory using the File-> Menu.




Go ahead and process the data by clicking on “Yes.” If you had previously worked with this directory and already processed it once, you can click “No.” This is useful when working with very large data sets you want to reanalyze.

Once the data is processed, click on the Ladder Assignment button and CAFA will automatically begin assigning the DNA ladder in your data:




Adjusting the sensitivity slider will adjust what CAFA considers to be a peak in the sequencing ladder. Normally, CAFA should be able to automatically guess all the peak positions correctly, but if the ladder trace is weak, going to a higher sensitivity can help.

CAFA defaults to the Beckman-400 ladder. If you are using a different ladder, you can edit the positions of the ladder peaks using the File-> Edit Ladder Positions dialog:




You can use different ladders for different traces, and each ladder is stored and associated with the particular trace you are looking at. If you do not have a ladder in your sample, CAFA may crash. The best is then to remove that data set from the directory, reprocess the data and continue the analysis.

Once you have assigned each ladder, CAFA will present you with the fitting parameters:




Here you can enter the left and right most peaks you would like to fit. If your Nucleic Acid is shorter than your ladder, you can select to fit between the peaks where there is data. Simply click on "Ok" to begin fitting.

With the data fit, CAFA will ask you where to save the raw data:



At this stage you have finished the fitting of your data, congratulations! It is usually not a good idea to save this file in your original data directory, as it is a txt file like the raw data and may interfere with the analysis if you run the program again on the raw data. Instead, creating a directory called Fit\_data separates the raw and fit data.




Above are the files that are now in your directory. The files starting with area\_F1 correspond the individual peak areas and have this format:

```
F1s1a_07071119A5.txt
60 159922.574379
61 83176.303827
62 78713.837047
63 109612.246825
64 59842.400391
65 49383.548448
66 80645.881063
67 75769.886078
68 66772.737142
69 64594.606567
70 57230.873266
71 29542.321073
72 50833.046553
73 44031.947376
74 28988.165851
75 13190.880051
76 195933.832520
```

where the first column is the numbering corresponding to the DNA ladder and the second column is the raw peak area. The Data file saved after all the fitting has a similar format with all the data concatenated into columns.

One final step can be carried out to normalize and filter the data. This step is only if you are using an indirect labeling scheme for your nucleic acid (e.g. reverse transcription). To do this select File-> Normalize and Filter. This can be done anytime after the fitting even after you have quit and restarted CAFA.

Select the Data file you saved after fitting and a Reference\_Selection window will popup:




Select from this window the name of your experiment containing the background trace (i.e. a trace where the RT reaction was run on a non-modified nucleic acid). In this case it is the area\_F1 seq1\_070711215M.dat data.

CAFA will then ask you where you want to save the normalized data, which you can save in the Fit\_data directory too. This data looks like this:

Residues	area_F1s1a_07071119A5.dat	area_F1s1b_070711215M.dat	area_F1seq1_070711215M.dat
60	3.356953e+00	4.242148e+00	1.542471e+00
61	1.745963e+00	2.387782e+00	1.208065e+00
62	1.652291e+00	2.154473e+00	7.477427e-01
63	2.300882e+00	2.748971e+00	1.886696e+00
64	1.256158e+00	1.471044e+00	3.605530e-01
65	1.036616e+00	1.372553e+00	3.705500e-01
66	1.692847e+00	1.984043e+00	1.586979e+00
67	1.590494e+00	1.938537e+00	6.243484e-01
68	1.401634e+00	1.817204e+00	4.868283e-01
69	1.355913e+00	1.574855e+00	1.571141e+00
70	1.201340e+00	1.370064e+00	1.309624e+00
71	6.201261e-01	6.345996e-01	0
72	1.067042e+00	1.309437e+00	9.398702e-01
73	9.242796e-01	1.109162e+00	2.474074e-01
74	6.084939e-01	6.001110e-01	0
75	NaN	NaN	NaN
76	NaN	NaN	NaN
77	NaN	NaN	NaN
78	3.094050e+00	3.699825e+00	1.870406e+00
79	NaN	NaN	NaN
80	NaN	NaN	NaN
81	NaN	NaN	NaN
82	5.732181e-01	6.727585e-01	1.625434e-01
83	5.701988e-01	6.333145e-01	3.072298e-01

Where the data is now normalized to the mean intensity and data near strong RT stops is labeled as NaN. This is non-reliable data since it is a result of an RT stop.





The above figure is plotted by CAFA to show where data was removed. CAFA also indicates the percentage of the data that was excluded (in this case 22.5%). You can adjust the sensitivity of the filtering by changing the Sensitivity (default 2) in the Reference\_Selection window.

At this stage, the data can be imported into excel and further analyzed.