



# Modeling and Simulation in Bulk

OpenSim Advanced Users Workshop August 16, 2011

# My Motivation



# Today's Example



#### You will need

- Matlab
- BatchOpenSimInMatlab.zip
  - Workshop Project Page

#### Tools in the Command Line

```
Microsoft Windows XP [Version 5.1.2600]
(C) Copyright 1985-2001 Microsoft Corp.

C:\Documents and Settings\emarnold\cd "My Documents\_ModelersKitchen\BulkProcessing\GIL06\IKSetup"

C:\Documents and Settings\emarnold\My Documents\_ModelersKitchen\BulkProcessing\GIL06\IKSetup\c:\OpenSim2.3.2\bin\ik.exe -S Setup_IK_GIL06_fast3.xml_
```

#### C:\OpenSim2.3.2\bin\ik.exe -S Setup\_File.xml

Path to OpenSim tool you want to use (ik.exe, analyze.exe, etc)

-S for setup

Setup File Name

## Output

```
_ | _ | ×
C:\WINDOWS\system32\cmd.exe
Microsoft Windows XP [Version 5.1.2600]
(C) Copyright 1985-2001 Microsoft Corp.
G:\Documents and Settings\emarnold>cd "My Documents\_ModelersKitchen\BulkProcess
ing\GIL06\IKSetup"
This program includes software developed by the
Apache Software Foundation (http://www.apache.org/).
Constructing tool from setup file Setup_IK_GIL06_fast3.xml.
Loaded model GIL06 from file ../GIL06_gait2392_simbody_v23.osim
            MODEL: GILØ6
           forces: 92
         analyses: 0
           bodies: 13
           .joints: 12
          markers: 39
Running tool GIL06_fast3.
Loaded marker file C:\Documents and Settings\emarnold\My Documents\ ModelersKitc
hen\BulkProcessing\GIL06\MarkerData\GIL06_fast3.trc (37 markers, 114 frames)
                      total squared error = 0.0223747, marker error: RMS=0.024
Frame 1 (t=0.125003):
5911. max=0.0602543 (toe_r)
                      total squared error = 0.0227277, marker error: RMS=0.024
Frame 2 (t=0.133337):
```

```
Frame 111 (t=1.04167): total squared error = 0.0299058, marker error: RMS=0.028
43, max=0.0641389 (thighwand_r)
Frame 112 (t=1.05): total squared error = 0.0221762, marker error: RMS=0.024
4818, max=0.0614205 (thighwand_r)
InverseKinematicsTool: 112 frames in 4.047s

IK compute time = 5063ms

C:\Documents and Settings\emarnold\My Documents\_ModelersKitchen\BulkProcessing\
GIL06\IKSetup>
```

captured in out.log

## Moving to Matlab

% write a string to send to the command line
Command = 'c:/OpenSim2.3.2/bin/ik.exe -S
 Setup\_File.xml';

% send it to the command line system(Command);

## Setup Files in Matlab


```
🖔 C:\Documents and Settings\emarnold\My Documents\_ModelersKitchen\BulkProcessing\GIL11\lKSetup\GIL11_gener... 🚐
File Edit Search View Format Language Settings Macro Run TextFX Plugins Window ?
  3 🖴 🗎 🖺 🖺 🧸 🨘 🚕 | 🚜 🖍 🖍 🏗 | Þ C | 🛗 🛬 | 🔍 🔇 🖫 🖫 🚍 🚍 🚛 🕦 🗐 🐷 🗩 🕪
GIL11_generic_Setup_IK_v23.xml
         <?xml version="1.0" encoding="UTF-8" ?>
       - < OpenSimDocument Version="20300">
             <InverseKinematicsTool name="GIL11 free2">
                 <!--Name of the .osim file used to construct a model.-->
                 <model file> ../GIL11 gait2392 simbody v23.osim </model file>
   6
                 <!--Markers and coordinates to be considered (tasks) and their weightings.-->
                 <IKTaskSet file="./gait2392 GIL IK Tasks.xml"/>
                 <!--TRC file (.trc) containing the time history of observations of marker
                     positions. -->
                 <marker file> .../MarkerData/GIL11 free1.trc </marker file>
  11
                 <!-- The name of the storage (.sto or .mot) file containing coordinate
  12
                     observations. Coordinate values from this file are included if there is
  13
                     a corresponding coordinate task .-->
  14
                 <coordinate file> Unassigned </coordinate file>
  15
                 <!--Time range over which the inverse kinematics problem is solved. -->
  16
                 <time range>
                                    0.13333 1.0917</time range>
  17
                 <!--Name of the motion file (.mot) to which the results should be written.-->
  18
                 <output motion file> ../IKResults/GIL11 free1 ik.mot </output motion file>
  19
                 <include markers>true</include markers>
                 <!--A positive scalar that is used to weight the importance of satisfying
  21
                     constraints. A weighting of 'Infinity' or if it is unassigned results
  22
                     in the constraints being strictly enforced. -->
  23
                 <constraint weight>
                                          20.00000000 </constraint weight>
  24
                 <!-- The accuracy of the solution in absolute terms. I.e. the number of
  25
                     significantdigits to which the solution can be trusted .-->
  26
                                  0.00001000 </accuracy>
  27
                 </InverseKinematicsTool>
  28
         </0penSimDocument>
  29
  30
eXtensible Markup nb char: 1681 nb line: 30
                                                 Ln:23 Col:66 Sel:0
                                                                                     Dos\Windows ANSI
```

## Setup Files in Matlab


- xmlread
  - Parse an xml file and turn it into a document object model
- Get and Set to interact with attributes
  - getElementsByTagName
  - getFirstChlid
  - setNodeValue
  - getAttributes
- xmlwrite
  - Write object back out to an xml file

#### Batch Example

- 2 Subjects
- 4 Gait Speeds
  - xslow, slow, free, fast
- Inverse Kinematics
- Muscle Analysis
- Plot results








https://simtk.org/home/mspeedwalksims

## Tips

- Give gait trials meaningful, consistent names during data collection.
  - ex. Run\_200\_02.trc not trial1.trc
- Avoid spaces in file and folder names
- Adopt regular folder structure for all subjects
- Use values and text from inputfiles rather than hard coding in matlab
- Verify input/output file conventions when moving between versions
- Use sendmail to automate updates

# Following up...



#### ...contribute to the Wiki!

